



Procurement Executive Ministry of Defence

CHEMICAL DEFENCE ESTABLISHMENT

Porton Down Salisbury Wilts SP4 OJQ

Telephone Idmiston 610211 ext 201
(STD 0980)

Inch

Prof M Meselson
Department of Biochemistry
and Molecular Biology
Harvard University
7 Divinity Avenue
Cambridge
Massachusetts 02138

Please reply to The Director
Your reference

Our reference

Date

7 September 1984

Dear Matt

Sorry for the delay in answering your letter of 6th August. Thanks for the information and samples. Analysis of the samples is in progress and I should have the results by the end of September.

I promised details of our procedures for work up and analysis. These I now enclose.

I note your November dates and will contact you about these later.

Yours sincerely

T D INCH

Enc

CDE ANALYSIS OF POLLEN AND ENVIRONMENTAL RESIDUES

Three methods are used depending on instrument availability, matrix and sensitivity required.


1. GC-MS (SIM) in EI mode using HFB derivs (VG 7070 EQ instrument).
2. GC-EC using HFB derivs, all trichothecenes detected as their hydrolysis products T-2 tetraol, scirpentriol, DON, NIV, 3,4,7,15-tetraol, verrucarol.
3. GC-MS (SIM) in -ve CI mode using PFP derivs (Finnegan 4500 quadrupole) [PFP derivs are used because MV & T-2 tetraol (HFB)₄ derivs give M⁻ ions beyond the 1000 amu mass range of the instrument. The VG 7070 EQ give rather poor sensitivity in -ve CI mode].

Methods 1 and 2 are frequently used in parallel. GC-EC was subject to too much interference when used to detect T-2 & DAS but for pollen matrixes it has been shown to be a reliable method for detecting the polyhydroxytrichothecenes. For other matrixes such as leaves GC-EC is subject to too much interference and GC-MS is used. Any positive by GC-EC would need to be confirmed by GC-MS. Best limits of detection are obtained with method 3 especially for polyhydroxytrichothecenes such as DON, NIV, T-2 tetraol, but is not routinely used because of non-availability of the instrument. Negative CI also has the advantages that it is less prone to interference from background. Method 3 is used mainly for blood and urine analysis where lack of detection ranging from 0.1-5 ppb can be obtained.

R M BLACK
CDE Porton Down
6/9/84

COE TRICHOCECENE ANALYSIS

POLLEN ANALYSIS - METHODS 1 & 2



METHOD 1

GC-MS CONDITIONS (HFB derivs, EI)

DEQUALIZATION - as for GC-EC

GC CONDITIONS

COLUMN : 25m BPI, 0.25mm dia

INJECTOR TEMP : 260°

OVEN TEMP : 160° for 1min, 160-275° at 10°/min, 275° for 5min

TRANSFER LINE : 260°

CARRIER : He, 15psi

SPLIT 1:1, 0.5 µl injected

MS CONDITIONS

INSTRUMENT : VG 7070 EQ

MODE : SIR, low resolution, EI

Two groups of 10 ions monitored, A 5 min 15 sec → 9 min 40 sec
B 9 min 48 sec → 16 min 30 sec

IONS MONITORED, RETENTION TIMES

Group	Ion	m/z	Retention Time	Group	Ion	m/z	Retention Time
<u>A</u>	NIV	1096, 1077	7 ³⁸	<u>B</u>	OAS	502, 674	11 ³⁵
	00N	884, 865	8 ³⁷		15-ALDON	730, 688	10 ²²
	sc-trid	870, 855	8 ²⁷		HT-2	655, 672	12 ⁵⁹
	T-2 tetracl	869, 868	8 ⁴⁵		T-2	519, 501	15 ²⁸
	FUS-X	942, 923	9 ³³		15-MAS	656, 625	9 ⁵⁶

Higher resolution may help for extracts of leaves where considerable background occurs.

Detection limit usually around 0.5ng - 5ng per sample depending on the background.

GC-EC CONDITIONS

DEGRATIZATION

1. To dry residue in 1 ml reaction vial add 200 µl toluene / oxygen (95:5)
2. Add 20 µl HPAI Cuddihill
3. Vortex 30 sec
4. Heat at 60° , 2h , cool.
5. Add 500 µl 5% NaHCO₃
6. Vortex 2 min.
7. Remove aqueous layer
8. Vortex 30 sec with 500 µl water (x 2) , remove aqueous layer.
9. Inject 0.5 µl of toluene soln.

GC CONDITIONS


INSTRUMENT : PYE 304
 COLUMN : 25m BPL , 0.25 dia.
 INJECTOR TEMP : 260°
 DETECTOR TEMP : 275°
 OVEN TEMP : 170° 5min , 170-190° 5°/min ; 190° 5min ;
 190-250° 4°/min ; 250° 5min.
 CARRIER : H₂ , 9 psi
 MAKE-UP : N₂ , 10 psi
 CATAL (Flow 20 ml/min)
 SPLIT FLOW : 10 ml/min
 SEPTUM PURGE : 3 ml/min

RETENTION TIME :	T-2 tetracl	14.8 min
scirpentriol	13.2 min	
nivalenol	11.0 min	
DNV	13.9 min	
3,4,7,15-tetracl	10.7 min	
verrucareol	14.6 min	

Occasional interfering peaks are usually separated using isothermal conditions.

Approx. detection limits : T-2, OAS, NEOS 5 ng per 1-10 mg pollen
 VER-A, VER-B, MV, DON 20-25 ng per 10 mg "

TRICHOHECENES: DETECTION IN BLOOD, URINE & POLLEN



concentrate

C₁₈SEP-PAK

H₂O eluant →

40% MeOH eluant
DON, NIV, Sc-triol
T-2 tetraol

80% MeOH eluant
T-2, HT-2, DAS
NEOS, Ac-DON

N.B. Toxins are fractionated so that the esterified trichothecenes can be extracted, thus avoiding any partial hydrolysis on concentration.

CH₂Cl₂ extract

Dry residue

Derivatize as PFP esters [PFPI 60°, 1h]

GC-MS, SIM, NCI (CH₄ reagent gas)

[detection limits for pollen ca. 0.1 - 5 ng per sample]